

# AC-81MD

AC Adaptor for Medical Devices



## Overview

### **Provides DC power for compatible Sony medical equipment**

The AC-81MD is an AC adaptor that is intended for use with these compatible Sony medical products\* only:

MCC-3000MT, NU-IP40S, NU-IP40D, UP-D711MD (AC-81MD available as optional accessory)

HVO-500MD/500MD(Surgical Version), HVO-550MD/550MD(Surgical Version), HMS-3000MT (AC-81MD included as standard)

The adaptor can operate wherever an AC mains outlet (100-240V, 50/60Hz) is available. It provides a regulated 24V DC output (3.3A max.) suitable for powering the products as detailed above.

Includes AC power cable and DC cable (1.8m) for connection to device.

\*Product availability may vary according to the sales area. For more information, please contact your nearest Sony sales office or an authorised dealer.

## Specifications

Specifications	
Power	AC IN:100 V - 240 V、 50/60 Hz, 1.0 A - 0.5 A DC OUT: 24 V 3.3 A
Operating Temperature	0 °C to 40 °C (32°F to 104 °F)
Operating Humidity	20 % ~ 85 % (no condensation allowed)
Storage/Transport Temperature	-20 °C to +60 °C (-4 °F to +140 °F)
Storage/Transport Humidity	20% to 90% (no condensation allowed)
Operating/Storage/Transport Pressure	700 hPa ~ 1,060 hPa
Dimensions (W x H x D) *1	60 × 35 × 160 mm (Not incl. projecting parts) 2 3/8 × 1 7/16 × 6 3/8 inches (Not incl. projecting parts)

Mass	400 g (AC adaptor only) 14 oz. (AC adaptor only)
------	---

---

Supplied Accessories	Instructions for Use (1) Service Contact list (1) AC power code (supplied at local) (1)
----------------------	---

---

## Notes

*1	The values for dimensions are approximate.
----	--

---

## Gallery

